

Anti-CLEC10A antibody (60-180) (STJ113856)

STJ113856

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10462
Gene Symbol	CLEC10A
Uniprot ID	CLC10_HUMAN
Immunogen	
Immunogen Region	60-180
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 60-180 of human CD301/CLEC10A (NP_878910.1).
Immunogen Sequence	FQNSKFQRDLVTLRTDFSNF TSNTVAEIQALTSQGSLEE TIASLKAIEVEGFKQERQAGV SELQEHTTQKAHLGHCPHCP SVCVPVHSEMLLRVQQLVQD LKKLTCQVATLNNNASTEGT C



Western blot analysis of various lysates using CD301/CLEC10A antibody (STJ113856) at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/ proteins: 25 Mu g per lane. Blocking buffer: 3 % nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081