

Anti-ABCD1 antibody (340-500) (STJ114034)

STJ114034

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/ICC:1:50-1:200
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	215
Gene Symbol	ABCD1
Uniprot ID	ABCD1_HUMAN
Immunogen	
Immunogen Region	340-500
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 340-500 of human ABCD1 (NP_000024.2).
Immunogen Sequence	SASGLLMVAVPIITATGYSE SDAEAVKKAALKKEEELVS ERTEAFTIARNLLTAAADAI ERIMSSYKEVTELAGYTARV HEMFQVFEDVQRCHFKRPRE LEDAQAGSGTIGRSGVRVEG PLKIRGQVVDVEQGIICENI PIVTPSGEWWASLNIRVEE G



Immunofluorescence analysis of NIH/3T3 cells using ABCD1 Rabbit polyclonal antibody (STJ114034) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of various lysates using ABCD1 Rabbit polyclonal antibody (STJ114034) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.