

Anti-SNRPA1 antibody (1-220) (STJ114054)
STJ114054

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB: 1:500-1:2000 IHC-P: 1:50-1:200 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6627
Gene Symbol	SNRPA1
Uniprot ID	RU2A_HUMAN
Immunogen	
Immunogen Region	1-220
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human SNRPA1 (NP_003081.2).
Immunogen Sequence	MVKLTAEIEQAQYTNVAVR DRELDLRGYKIPVIENLGAT LDQFDAIDFSDNEIRKLDGF PLLRRLKTLVNNNRICRIG EGLDQALPCLTELITNNSL VELGDLPLASLKSITYLSI LRNPVTNKKHYRLYVIYKVP QVRVLDQKVKLKERQEAEK MFKGKRGQAQLAKDIARRSKT FNPAGALPTDKKKGPPSPGD VEAIKNAIANASTLAEVERL




Immunohistochemistry analysis of SNRPA1 in paraffin-embedded rat lung using SNRPA1 Rabbit polyclonal antibody (STJ114054) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of SNRPA1 in paraffin-embedded rat liver using SNRPA1 Rabbit polyclonal antibody (STJ114054) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of SNRPA1 in paraffin-embedded mouse liver using SNRPA1 Rabbit polyclonal antibody (STJ114054) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of SNRPA1 in paraffin-embedded human lung cancer using SNRPA1 Rabbit polyclonal antibody (STJ114054) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081