

Anti-AGMAT antibody (1-352) (STJ114084)

STJ114084

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79814
Gene Symbol	AGMAT
Uniprot ID	GDAH_HUMAN
Immunogen	
Immunogen Region	1-352
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-352 of human AGMAT (NP_079034.3).
Immunogen Sequence	MLRLLASGCARGPGPGVGAR PAAGLFHPGRRQSRQASDAP RNQPPSPFVARPVGVCSMM RLPVQTSPEGLDAAFIGVPL DTGTSNRPGARFGPRRIREE SVMLGTVNPSTGALPFQSLM VADLGDVNVNLYNLQDSCRR IQEAYEKIVAAGCIPLTLGG DHTITYPILQAMAKKHGPGV LLHVDAHTDITDKALGEKLY HGAPFRRCVDEGLLDCKRWV QIGIRGSSTTLDPYRYNRS



Western blot analysis of various lysates, using AGMAT Rabbit polyclonal antibody (STJ114084) at 1:2000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 45s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081