

Anti-PCSK1 antibody (28-170) (STJ114128)

STJ114128

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5122
Gene Symbol	PCSK1
Uniprot ID	NEC1_HUMAN
Immunogen	
Immunogen Region	28-170
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 28-170 of human PCSK1 (NP_000430.3).
Immunogen Sequence	KRQFVNEWAAEIPGGPEAAS AIAEELGYDLLGQIGSLENH YLFKHKNHPRRSRRSAFHIT KRLSDDDRVIWAEQQYEKER SKRSALRDSALNLFNDPMWN QQWYLDTRMTAALPKLDLH VIPVWQKGITGKGWITVLD DGL



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081