

Anti-PPARGC1A antibody (610-798) (STJ114227)

STJ114227

GENERAL INFORMATION





Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	10891
Gene Symbol	PPARGC1A
Uniprot ID	PRGC1_HUMAN
Immunogen	
Region	610-798
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 610-798 of human PGC1 Alpha (NP_037393.1).
Immunogen Sequence	HRTHRNSPLYVRSRSPYS RRPYDSYEEYQHERLKREE YRREYKRESEKQREKQK QKAIERRVIYVGKIRPDTT RTELDRFEVFGIEIECTVN LRDDGDSYGFTYRYTCDAF AALENGYTLRRSNETDFELY FCGRKQFFKSNYADLDSNSD DFDPASTKSKYDSDLDFDSSL KEAQRSLRR



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081