

Anti-CALU antibody (20-315) (STJ114282)

STJ114282

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:100
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	813
Gene Symbol	CALU
Uniprot ID	CALU_HUMAN
Immunogen	
Immunogen Region	20-315
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-315 of human CALU (NP_001210.1). KPTKKDRVHHEPQLSDKVH NDAQSFDYDHDAFLGAEEAK TFDQLTPEESKERLGKIVSK IDGDKDGFVTVDELKDWIKF AQKRWIYEDVERQWKGHDLN EDGLVSWEYKKNATYGYVLD DPDPDDGFNYKQMMVRDERR FKMADKDGDLIATKEEFTAF LHPPEYDYMKDIVVQETMED IDKNADGFIDLEEYIGDMYS HDGNTDEPEWVKTEREQFVE FRDKNRDGMKDKEETKDWI
Immunogen Sequence	



Immunofluorescence analysis of HeLa cells using CALU antibody (STJ114282). Blue: DAPI for nuclear staining.