

Anti-ARHGEF2 antibody (1-350) (STJ114421)

STJ114421

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9181
Gene Symbol	ARHGEF2
Uniprot ID	ARHG2_HUMAN
Immunogen	
Immunogen Region	1-350
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human ARHGEF2 (NP_004714.2).
Immunogen Sequence	MKEAKDARYTNGHLFTTISV SGMTCYACNKSITAKEALI CPTCNVTIHNRCCKDTLANCT KVKQKQQKAALLKNNTALQS VSLRSKTTIRERPSSAIYPS DSFRQSLLGSRRGRSSLSLA KSVSTTNIAGHFNDESPGL RRILSQSTDLSNMRNRTL SV ESLIDEEVIYSELMSDFEMD EKDFAADSWSLAVDSSFLQQ HKKEVMKQQDVYELIQTEL HHVRTLKIMTRLFRTGMLE