

Anti-SAFB2 antibody (1-255) (STJ114424)

STJ114424

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	9667
Gene Symbol	SAFB2
Uniprot ID	SAFB2_HUMAN
Immunogen	
Immunogen Region	1-255
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-255 of human SAFB2 (NP_055464.1).
Immunogen Sequence	MAETLPGSGDSGPGTASLGP GVAETGTRRLSELRVLDLRA ELKKRNLDTGGNKSVLMERL KKAVKEEGQDPDEIGIELEA TSKKSARKRCVKGLKMEEEGT EDNGLEDDSRDGDQEDMEASL ENLQNMGMMDMSVLDTEVA NSSAPDFGEDGTDGLLDSFC DSKEYVAAQLRQLPAQPPEH AVDGEFGKNTLETSSLNFKV TPDIEESLLEPENEKILDIL GETCKSEPVKESSELEQP



Immunohistochemistry of paraffin-embedded human placenta using SAFB2 rabbit polyclonal antibody (STJ114424) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using SAFB2 rabbit polyclonal antibody (STJ114424) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse kidney using SAFB2 rabbit polyclonal antibody (STJ114424) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081