

Anti-MAD2L2 antibody (1-211) (STJ114433)
STJ114433

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10459
Gene Symbol	MAD2L2
Uniprot ID	MD2L2_HUMAN
Immunogen	
Immunogen Region	1-211
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-211 of human MAD2B/MAD2B/MAD2L2 (NP_006332.3).
Immunogen Sequence	MTTLTRQDLNFGQVADVLC EFLEVAVHLILYREVYPVG IFQKRKKYNNVPVQMSCHPEL NQYIQDTLHCVKPLLEKNDV EKWVVILDKEHRPVEKFVF EITQPPLLSISSDSLLSHVE QLLRAFILKISVCDVLDHN PPGCTFTVLVHTREAATRNM EKIQVIKDFPWILADEQDVH MHDPRLIPLKMTSDILKMQLYVEERAHKGS




Immunofluorescence analysis of U2OS cells using MAD2B/MAD2B/MAD2L2 antibody (STJ114433) at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of L929 cells using MAD2B/MAD2B/MAD2L2 antibody (STJ114433) at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of H9C2 cells using MAD2B/MAD2B/MAD2L2 antibody (STJ114433) at dilution of 1:100. Blue: DAPI for nuclear staining.



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081