

Anti-RASSF1 antibody (101-340) (STJ114441)

STJ114441

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11186
Gene Symbol	RASSF1
Uniprot ID	RASF1_HUMAN
Immunogen Region	101-340
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 101-340 of human RASSF1 (NP_009113.3).
Immunogen Sequence	CCGPRDLGWEPVERDTNVD EPVEWETPDLSQAEIEQKIK EYNAQINSNLFMSLNKDGSY TGFIVQLKLVVPVSPSSK KPPSLQDARRGPGRGTSVRR RTSFYLPKDAVKHLHVLSRT RAREVIEALLRKFLVDDPR KFALFERAERHGGVYLRKLL DDEQPLRLRLLAGPSDKALS FVLKENDSGEVNWDAFSMPE LHNFLRLQREEEHLRQIL QKYSYCRQKIQEALHACPL