

Anti-CHRFAM7A antibody (1-150) (STJ114477)

STJ114477

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	89832
Gene Symbol	CHRFAM7A
Uniprot ID	CRFM7_HUMAN
Immunogen	
Immunogen Region	1-150
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human CHRFAM7A (NP_647536.1). MQKYCIYQHFQFQLLIQHLW IAANCDIADERFDATFHTNV LVNSSGHCQYLPPGIFKSSC YIDVRWFPPFDVQHCKLKFGS
Immunogen Sequence	WSYGGWSLDLQMQEADISGY IPNGEWDLVGIPGKRSEFY ECCKEYPYDVTFTVTMRRT LYGLNLLIP



Western blot analysis of various lysates using CHRFAM7A Rabbit polyclonal antibody (STJ114477) at 1:400 dilution. Secondary antibody:HRP Goat Anti-Rabbit IgG (H+L) (STJS000858) at 1:10000 dilution. Lysates/ proteins: 25 µg per lane. Blocking buffer: 3 % nonfat dry milk in TBST. Detection:ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081