

Anti-TRPV3 antibody (400-500) (STJ114487)

STJ114487

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	162514
Gene Symbol	TRPV3
Uniprot ID	TRPV3_HUMAN
Immunogen	
Immunogen Region	400-500
Specificity	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human TRPV3 (NP_659505.1).
Immunogen Sequence	DNSVLEITVYNTNIDNRHEM LTLEPLHTLLHMKWKKFAKH MFFLSFCFYFFYNITLTVS YYRPREEEAIPHPLALTHKM GWLQLLGRMFVLIWAMCISV K