

Anti-HFM1 antibody (1-240) (STJ114488)

STJ114488

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	164045
Gene Symbol	HFM1
Uniprot ID	HFM1_HUMAN
Immunogen	
Immunogen Region	1-240
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human HFM1 (NP_001017975.4).
Immunogen Sequence	MLKSNDCLFSLLENLFFEKPD EVENHPDNEKSLDWFLPPAP LISEIPDTQEEEELESHKL LGQEKRPKMLTSNLKITNED TNYISLTQKFQFAFPSDKYE QDDLNLLEGVGNNDLSHVAGK LTYASQKYKNHIGTEIAPEK SVPDDTKLVNFAEDKGESTS VFRKRLFKISDNIHGSAYSN DNELDSHIGSVKIVQTEMNK GKSRNYSNSKQKFQYSANVF TANNAFSASEIGEGMFKAP



Western blot analysis of extracts of various cell lines, using HFM1 antibody (STJ114488) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 60s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081