

Anti-PAFAH1B1 antibody (1-170) (STJ114517)

STJ114517

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5048
Gene Symbol	PAFAH1B1
Uniprot ID	LIS1_HUMAN
Immunogen	
Immunogen Region	1-170
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human PAFAH1B1 (NP_000421.1). MVLSSQRQRDELNRAIADYLR SNGYEEAYSVFKKEAELDVN EELDKKYAGLLEKKWTSVIR LQKKVMELESKLNEAKEEFT
Immunogen Sequence	SGGPLGQKRDPEWIPRPPE KYALSGHRSPVTRVIFHPVF SVMVSASEDATIKVWDYETG DFERTLKGHTDSVQDISFDH SGKLLASCSA



Western blot analysis of extracts of various cell lines, using PAFAH1B1 antibody (STJ114517) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081