

Anti-NUP35 antibody (1-326) (STJ114635)

STJ114635

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	129401
Gene Symbol	NUP35
Uniprot ID	NUP35_HUMAN
Immunogen	
Immunogen Region	1-326
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-326 of human NUP35 (NP_612142.2).
Immunogen Sequence	MAAFAVEPQGPALGSEPMMML GSPTSPKPGVNAQFLPGFLM GDLPAVTPQPRSIGPSVG VMEMRSPLLAGGSPQPVPV AHKDKSGAPPVRSIYDDISS PGLGSTPLTSRRQPNISVMQ SPLVGVSTPTGTGQSMFSPA SIGQPRKTTLSPAQLDPFYT QGDSLTSDEHLDDSWVTVFG FPQASASYILLQFAQYGNIL KHVMSNTGNWMHIRYQSKLQ ARKALSKDGRIFGESIMIG



Western blot analysis of various lysates using NUP35 Rabbit polyclonal antibody (STJ114635) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081