

Anti-FEM1B antibody (368-627) (STJ114668)

STJ114668

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	10116
Gene Symbol	FEM1B
Uniprot ID	FEM1B_HUMAN
Immunogen	
Immunogen Region	368-627
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 368-627 of human FEM1B (NP_056137.1).
Immunogen Sequence	LHALHLRQKGNRNTHKDLLR FAQVFSQMIHLNETVKAPDI ECVLRCSVLEIEQSMNRVKVNS ISDADVHNAMDNIECNLYTF LYLVCISTKTQCSEEDQCKI NKQIYNLIHLDPRTREGFTL LHLAVNSNTPVDDFHTNDVC SFPNALVTKLLDCCGAEVNA VDNEGNSALHIIVQYNRPIS DFLTLHSIIISLVEAGAHTD MTNKQNKTPDKSTTGVSSEI LLKTQMKMSLKCLAAARAVR



Western blot analysis of extracts of A-549 cells, using FEM1B antibody (STJ114668) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081