

Anti-NKAIN4 antibody (89-208) (STJ114707)

STJ114707

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	128414
Gene Symbol	NKAIN4
Uniprot ID	NKAI4_HUMAN
Immunogen	
Immunogen Region	89-208
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 89-208 of human NKAIN4 (NP_690603.3).
Immunogen Sequence	LLKDSSELLTFSLSRHRSWWERWPGCLHEEVPAVGLGAPH GQALVSGAGCALEPSYVEAL HSCLQILIALLGFCGCGQVV SVFTEEDSFDFIGGFDFPPF LYHVNEKPSSLLSKQVYLPA