

# Anti-SLC27A1 antibody (510-646) (STJ114713)

STJ114713

## GENERAL INFORMATION


|                          |                    |
|--------------------------|--------------------|
| <b>Product Type</b>      | Primary antibodies |
| <b>Short Description</b> |                    |
| <b>Applications</b>      | WB/IF/ICC/ELISA    |
| <b>Host/Source</b>       | Rabbit             |
| <b>Reactivity</b>        | Human/Mouse/Rat    |

## PRODUCT PROPERTIES


|                            |                                                                                                                                        |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                                                             |
| <b>Clone ID</b>            |                                                                                                                                        |
| <b>Concentration</b>       | Lot specific                                                                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                                                           |
| <b>Purification</b>        | Affinity purification                                                                                                                  |
| <b>Dilution</b>            | WB:1:1000-1:2000                                                                                                                       |
| <b>Range</b>               | IF/CC:1:50-1:200                                                                                                                       |
|                            | ELISA:Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Formulation</b>         | PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.                                                                                     |
| <b>Isotype</b>             | IgG                                                                                                                                    |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                         |

## TARGET INFORMATION


|                           |                                                                                                                                               |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 376497                                                                                                                                        |
| <b>Gene Symbol</b>        | SLC27A1                                                                                                                                       |
| <b>Uniprot ID</b>         | S27A1_HUMAN                                                                                                                                   |
| <b>Immunogen</b>          |                                                                                                                                               |
| <b>Immunogen Region</b>   | 510-646                                                                                                                                       |
| <b>Specificity</b>        | Recombinant fusion protein containing a sequence corresponding to amino acids 510-646 of human SLC27A1 (NP_940982.1).                         |
| <b>Immunogen Sequence</b> | DTFRWRGENVSTTEVEGVLS RLLGQTDVAVYGVAVPGVEG KAGMAAVADPHSLDPNAYI QELQKVLAPYARPIFLRLLP QVDTTGTFKIQKTRLQREGF DPRQTSDDLFLDLKQGHYL PLNEAVYTRICSGAFAL |




Immunofluorescence analysis of U2OS cells using SLC27A1 Rabbit polyclonal antibody (STJ114713) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of PC-12 cells using SLC27A1 Rabbit polyclonal antibody (STJ114713) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH/3T3 cells using SLC27A1 Rabbit polyclonal antibody (STJ114713) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of HeLa cells using SLC27A1 Rabbit polyclonal antibody (STJ114713) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081