

**Anti-LYPD1 antibody (21-141) (STJ114740)**

STJ114740

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB:1:1000-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">116372</a>
<b>Gene Symbol</b>	<a href="#">LYPD1</a>
<b>Uniprot ID</b>	<a href="#">LYPD1_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	21-141
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 21-141 of human LYPD1 (NP_653187.3).
<b>Immunogen Sequence</b>	LQIQCYQC EEFQLNND CSSP EFIVNCTVNVQDMCQKEVME QSAGIMYRKSCASSAACLIA SAGYQSF CSPGKLNSVCISC CNTPLCNGPRPKKRGSSASA LRPGLRTTILFLKLALFSAH C