

Anti-KCTD12 antibody (196-325) (STJ114783)

STJ114783

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	115207
Gene Symbol	KCTD12
Uniprot ID	KCD12_HUMAN
Immunogen	
Immunogen Region	196-325
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 196-325 of human KCTD12 (NP_612453.1).
Immunogen Sequence	TPSQLDGSRRSGYITIGYR GSYTIGRDAQADAKFRRVAR ITVCGKTS LAKEVFGDTLNE SRDPDRPPERYTSRYYLKFN FLEQAFDKLSESGFHMVACS STGTCAFASSTDQSEDKIWT SYTEYVFCRE