

Anti-FUT8 antibody (376-575) (STJ114799)

STJ114799

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2530
Gene Symbol	FUT8
Uniprot ID	FUT8_HUMAN
Immunogen	
Immunogen Region	376-575
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 376-575 of human FUT8 (NP_835367.1). FHPIEEYMVHVEEHFQLLAR RMQVDKKRVYLATDDPSLLK EAKTKYPNYEFISDNSISWS AGLHNRYTENSRLRGVILDIH FLSQADFLVCTFSSQVCRVA YEIMQTLHPDASANFSLDD IYYFGGQNAHNQIAIYAHQP RTADEIPMEPGDIIGVAGNH WDGYSGVNRKLGRTGLYPS YKVREKIETVKYPTYPEAEK
Immunogen Sequence	



Western blot analysis of various lysates using FUT8 Rabbit polyclonal antibody (STJ114799) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081