

Anti-SEPT5 antibody (1-369) (STJ114819)

STJ114819

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	IHC-P:1:100-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	5413
Gene Symbol	SEPTIN5
Uniprot ID	SEPT5_HUMAN
Immunogen	
Immunogen Region	1-369
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-369 of human Septin 5 (NP_002679.2).
Immunogen Sequence	MSTGLRYKSKLATPEDKQDI DKQYVGFATLPNQVHRKSVK KGFDFTLMVAGESGLGKSTL VHSLFLTDLYKDRKLLSAEE RISQTVELKHTVDIEEKG V KLKLTIVDTPGFGDAVNNT CWKPTIDYVDQQFEQYFRDE SGLNRKNIQDNRVHCCLYFI SPFGHGLRPVDVGFMKALHE KVNIVPLIAKADCLVPSEIR KLKERIREEDKFGIHVYQF PECDSDEDEDFFKQDDRELK




Immunofluorescence analysis of U-2 OS cells using Septin 5 Rabbit polyclonal antibody (STJ114819) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunohistochemistry analysis of Septin 5 in paraffin-embedded mouse brain using Septin 5 Rabbit polyclonal antibody (STJ114819) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of Septin 5 in paraffin-embedded rat brain using Septin 5 Rabbit polyclonal antibody (STJ114819) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Western blot analysis of various lysates, using Septin 5 Rabbit polyclonal antibody (STJ114819) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081