

Anti-EBF1 antibody (500-591) (STJ114934)

STJ114934

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1879
Gene Symbol	EBF1
Uniprot ID	COE1_HUMAN
Immunogen	
Immunogen Region	500-591
Specificity	A synthetic peptide corresponding to a sequence within amino acids 500-591 of human EBF1 (NP_076870.1).
Immunogen Sequence	GGSPITFLNGSAANSPIAIVP SSPTMASSTSLPSNCSSTSG IFSFSPANMVS AVKQKSAFA PVVRPQTSPPTCTSTNGNS LQAISGMIVPPM



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081