

Anti-VPS39 antibody (606-875) (STJ115049)

STJ115049

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:100-1:500
Range	IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23339
Gene Symbol	VPS39
Uniprot ID	VPS39_HUMAN
Immunogen	
Immunogen Region	606-875
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 606-875 of human VPS39 (NP_056104.2).
Immunogen Sequence	RFHNCILQLYCEKVQGLMKE YLLSFPAGKTPVPAGEEEGE LGEYRQKLLMFLEISSYYDP GRLICDFPFDGLLEERALLL GRMGKHEQALFIYVHILKDT RMAEEYCHKHYDRNKDGNKD VYLSLLRMYLSPPSIHCLGP IKLELLEPKANLQAALQVLE LHHSKLDTTKALNLLPANTQ INDIRIFLEKVLLENAQKKR FNQVLKNLLHAEFLRVQEER ILHQQVKCIITEEKVCMVC



Western blot analysis of lysates from Rat brain, using VPS39 Rabbit polyclonal antibody (STJ115049) at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000896) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081