

Anti-AADAT antibody (34-429) (STJ115056)

STJ115056

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	51166
Gene Symbol	AADAT
Uniprot ID	AADAT_HUMAN
Immunogen	
Immunogen Region	34-429
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 34-429 of human AADAT (NP_001273611.1).
Immunogen Sequence	PKSMISLAGGLPNPNMFPEK TAVITVENGKTIQFGEEMMK RALQYSPSAGIPELLSWLKQ LQIKLHNPPTIHPPSQGQM DLCVTSQSGQQGLCKVFEMII NPGDNLVLLDEPAYSGTLQSL HPLGCNIIIVASDESGIVPD SLRDILSRWKPEDAKNPQKN TPKFLYTVPNGNNPTGNSLT SERKKEIYELARKYDFLIIE DDPYYFLQFNKFRVPTFLSM DVDGRVIRADSFSKIISG



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081