

Anti-DCAF13 antibody (298-597) (STJ115179)

STJ115179

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	25879
Gene Symbol	DCAF13
Uniprot ID	DCA13_HUMAN
Immunogen	
Immunogen Region	298-597
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 298-597 of human DCAF13 (NP_056235.4).
Immunogen Sequence	EEEP LHTILGKTVYTGIDHH WKEAVFATCGQQVDIWDEQR TNPICSM TWGFDSISSVKFN PIETFLLGSCASDRNIVLYD MRQATPLKKVILDMRTNTIC WNPMEAFIFTAANEDYNLYT FDMRALDTPVMVHMDHVS AV LDVDYSPTGKEFVSASF DKS IRIFPVDKSRSREVYHTKRM QHVICVKWTSDSKYIMCGSD EMNIRLWKANASEKLGVLTS REKAAKDYNQKLKEKFQHY



Western blot analysis of various lysates using DCAF13 Rabbit polyclonal antibody (STJ115179) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081