

Anti-DSP antibody (2609-2830) (STJ115263)
STJ115263

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1832
Gene Symbol	DSP
Uniprot ID	DESP_HUMAN
Immunogen	
Immunogen Region	2609-2830
Specificity	A synthetic peptide corresponding to a sequence within amino acids 2609-2830 of human Desmoplakin (NP_004406.2).
Immunogen Sequence	FSDTLEESSPIAIFDTENL EKISITEGIERGIVDSITGQ RLLEAQAQTGGIHPPTTGQK LSLQDAVSQGVVDQDMATRL KPAQKAFIGFEGVKGKKKMS AAEAVKEKWLPYEAGQRFLE FQYLTGGLVDPEVHGRISTE EAIRKGFIDGRAAQRLQDTS SYAKILTCPKTKLISYKDA INRSMVEDITGLRLLEASV SSKGLPSPYNMSSAPGSRSG SR



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081