

Anti-FLT4 antibody (25-776) (STJ115268)

STJ115268

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2324
Gene Symbol	FLT4
Uniprot ID	VGFR3_HUMAN
Immunogen	
Immunogen Region	25-776
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1174-1363 of human VEGFR3/FLT4 (NP_891555.2).
Immunogen Sequence	QGRGLQEEEEVCMAPRSSQS SEEGSFSQVSTMALHIAQAD AEDSPPSLQRHSLAARYYNW VSFPGCLARGAETRGSRRMK TFEEFPMPTPTYKGSVDNQT DSGMVLASEEFQIESRHRQ ESGFSCCKGPGQNVAVTRAHP DSQGRRRRPERGARGGQVFY NSEYGELSEPSEEDHCPSA RVTFFTDNSY



Immunofluorescence analysis of MCF7 cells using VEGFR3/FLT4 Rabbit polyclonal antibody (STJ115268) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081