

Anti-IL12RB1 antibody (25-305) (STJ115283)

STJ115283

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3594
Gene Symbol	IL12RB1
Uniprot ID	I12R1_HUMAN
Immunogen	
Immunogen Region	25-305
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 25-305 of human IL12RB1 (NP_714912.1).
Immunogen Sequence	RTSECCFQDPPYPDADSGSA SGPRDLRCYRISSDRYECSW QYEGPTAGVSHFLRCCLSSG RCCYFAAGSATRLQFSDQAG VSVLYTVTLWVESWARNQTE KSPEVTLQLYNSVKYEPLG DIKVSCLAGQLRMEWETPDN QVGAEVQFRHRTSPSPWKLG DCGPQDDDTESCLCPLEMNV AQEFQLRRRQLGSQGSSWSK WSSPVCVPPENPPQPQVRFS VEQLGQDGRRLTLKEQPT