

Anti-PDK4 antibody (102-411) (STJ115300)

STJ115300

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	5166
Gene Symbol	PDK4
Uniprot ID	PDK4_HUMAN
Immunogen	
Immunogen Region	102-411
Specificity	A synthetic peptide corresponding to a sequence within amino acids 102-411 of human PDK4 (NP_002603.1).
Immunogen Sequence	FHEKSPDDQKALSDFVDTLI KVRNRHHNVPTMAQGIIEY KDACTVDPVTNQNLQYFLDR FYMNRISTRMLMNQHILIFS DSQTGNPSHIGSIDPNCDDV AVQDAFECSRMLCDQYYLS SPELKTQVNGKFPDQPIHI VVPSHLHMLFELFKNAMR ATVEHQENQPSLTPIEVIV LGKEDLTIKISDRGGGVPLR IIDRLFSYTYSTAPTPVMDN SRNAPLAGFGYGLPISRLY




Immunofluorescence analysis of U2OS cells using PDK4 Rabbit polyclonal antibody (STJ115300) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunohistochemistry analysis of paraffin-embedded mouse heart using PDK4 Rabbit polyclonal antibody (STJ115300) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human placenta using PDK4 Rabbit polyclonal antibody (STJ115300) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded rat heart using PDK4 Rabbit polyclonal antibody (STJ115300) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081