

Anti-DNAJA3 antibody (221-480) (STJ115348)

STJ115348

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9093
Gene Symbol	DNAJA3
Uniprot ID	DNAJA3_HUMAN
Immunogen	
Immunogen Region	221-480
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 221-480 of human DNAJA3 (NP_005138.3).
Immunogen Sequence	AKGVNKEFTVNIMDTCCRCN GKGNEPGTKVQHCHYCGGSG METINTGPFVVRSTCRRCGG RGSIIISPCVVCRGAGQAKQ KKRVMIPVPAGVEDGGQTVRM PVGKREIFITFRVQKSPVFR RDGADIHSDFISIAQALLG GTARAQGLYETINVTIPPGT QTDQKIRMGGKGIPRINSYG YGDHYIHIKIRVPKRLTSRQ QSLILSYAEDETDVEGTVNG VTLTSSGGSTMDSSAGSKA



Immunofluorescence analysis of U2OS cells using DNAJA3 antibody (STJ115348). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of C6 cells using DNAJA3 antibody (STJ115348) at dilution of 1:100. Blue: DAPI for nuclear staining.