

Anti-CFHR3 antibody (80-200) (STJ115362)

STJ115362

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10878
Gene Symbol	CFHR3
Uniprot ID	FHR3_HUMAN
Immunogen	
Immunogen Region	80-200
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 80-200 of human CFHR3 (NP_066303.2).
Immunogen Sequence	AVPCLRKCYFPYLENGYNQN YGRKFVQGNSTEVACHPGYG LPKAQTTVTCTEKGWSPTPR CIRVRTCSKSDIEIENGFI
Sequence	ESSSIYILNKEIQYKCKPGY ATADGNSSGSITCLQNGWSA Q