

Anti-ZNF346 antibody (1-294) (STJ115368)
STJ115368

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	23567
Gene Symbol	ZNF346
Uniprot ID	ZN346_HUMAN
Immunogen	
Immunogen Region	1-294
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-294 of human ZNF346 (NP_036411.1).
Immunogen Sequence	MEYPAPATVQAADGGAAGPY SSSELLEGQEPDGVRFDRER ARRLWEAVSGAQPVGREEVE HMIQKNQCLFTNTQCKVCCA LLISESQKLAHYQSKKHANK VKRYLAIHGMETLKGTEKKL DSDQKSSRSKDNQCCPICN MTFSSPVVAQSHYLGKTHAK NLKLLKQQSTKVEALHQNREM IDPDKFCSLCHATFNDPVMA QQHYVGGKKHRKQETKLLMA RYGRADPAVTDFFPAGKGY



Immunohistochemistry of paraffin-embedded rat brain using ZNF346 antibody (STJ115368) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human tonsil using ZNF346 antibody (STJ115368) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using ZNF346 antibody (STJ115368) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081