

Anti-VPS4A antibody (1-120) (STJ115377)

STJ115377

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	27183
Gene Symbol	VPS4A
Uniprot ID	VPS4A_HUMAN
Immunogen	
Immunogen Region	1-120
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human VPS4A (NP_037377.1).
Immunogen Sequence	MTTSTLQKAIDLVTKATEED KAKNYEEALRLYQHAVEYFL HAIKYEAHSDKAKESIRAKC VQYLDRAEKLKDYLRSEKEH GKKPVKENQSEGKGSDDSE GDNPEKKKLQEQLMGAVVME