

Anti-ANGPTL4 antibody (26-406) (STJ115386)

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51129
Gene Symbol	ANGPTL4
Uniprot ID	ANGL4_HUMAN
Immunogen	
Region	26-406
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 26-406 of human ANGPTL4 (NP_647475.1).
Immunogen Sequence	GPVQSKSPRFASWDEMNVL A HGLLQLGQGLREHAERTRSQ LSALERRLSACGSACQGTEG STDPLAPESRVDPEVLHSL QTQLKAQNSRIQQLFHKVAQ QQRHLEKQHLRIQHLQSQFG LLDHKHLDHEVAKPARRKRL PEMAQPVDPAHNVSRHLRLP RDCQELFQVGERQSGLFELQ PQGSPPFLVNCKMTSDGGWT VIQRRHDGSVDFNRPWEAYK AGFGDPHGEFWLGLEKVHS