

Anti-DIABLO antibody (1-232) (STJ115393)

STJ115393

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	56616
Gene Symbol	DIABLO
Uniprot ID	DBLOH_HUMAN
Immunogen	
Immunogen Region	1-232
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-232 of human Smac/Smac/Diablo (NP_063940.1).
Immunogen Sequence	MAALKSWLSRSVTSFFRYRQ CLCVPVVANFKKRCFSELIR PWHKTVTIGFVTLCAVPIA QKSEPHSLSSSEALMRRRAVSL VTDSTSTFLSQTTYALIEAI TEYTKAVYTLTSLYRQYTS L LGKMNSEEEDEWVQVIIGAR AEMTSKHQEYLKLETTWMTA VGLSEMAAEAAAYQTGADQAS ITARNHIQLVKLQVEEVHQL SRKAETKLAEAEELRQKT QEEGEERAASEQ



Immunofluorescence analysis of MCF-7 cells using Smac/Smac/Diablo antibody (STJ115393). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081