

Anti-NSUN6 antibody (220-469) (STJ115415)
STJ115415

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	221078
Gene Symbol	NSUN6
Uniprot ID	NSUN6_HUMAN
Immunogen	
Immunogen Region	220-469
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 220-469 of human NSUN6 (NP_872349.1).
Immunogen Sequence	NLP ₃ SLVSHVLPQPGKIL DLCAAPGGKTHIAALMHDQ GEVIALDKIFNKVEIKQNA LLLGLNSIRAFCDGKAVK LDMVEDTEGEPFLPESFDR ILLDAPCSGMGQRPNMACTW SVKEVASYPQLQRLFTA ₃ AV QLLKPEGVLVYSTCTITLAE NEEQVAWALT ₃ KFPCLQLQPQ EPQIGGEGMRGAGLSCEQLK QLQRFDP ₃ SAVPLPDTMDSL REARREDMLRLANKDSIGF



Immunofluorescence analysis of MCF-7 cells using NSUN6 antibody (STJ115415). Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using NSUN6 antibody (STJ115415) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ3000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081