

Anti-IL7R antibody (21-239) (STJ115463)

STJ115463

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3575
Gene Symbol	IL7R
Uniprot ID	IL7RA_HUMAN
Immunogen	
Immunogen Region	21-239
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 21-239 of human CD127/IL7R (NP_002176.2).
Immunogen Sequence	ESGYAQNGDLEDAELDDYSF SCYSQLEVNQSQHSLTCAFE DPDVNITNLEFEICGALVEV KCLNFRKLQEIFYETKKFL LIGSNICVKVGEKSLTCKK IDLTITVKPEAPFDLSVVYR EGANDFVVTFTNTSHLQKKYV KVLMDVAYRQEKDENKWTW VNLSTKLTLLQRKLQPAAM YEIKVRSIPDHYFKGFWESEW SPSYYFRTPEINNSSGEMD



Western blot analysis of extracts of mouse spleen, using CD127/IL7R antibody (STJ115463) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081