

Anti-RPL11 antibody (1-177) (STJ115503)

STJ115503

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6135
Gene Symbol	RPL11
Uniprot ID	RL11_HUMAN
Immunogen	
Immunogen Region	1-177
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-177 of human RPL11 (NP_001186731.1).
Immunogen Sequence	MADQGEKENPMRELRIKLC LNICVGESGDLTRAQVLE QLTGQTPVFSKARYTVRSFG IRRNEKIAVHCTVRGAKAEE ILEKGLKVREYELRKNFSD TGNFGFGIQEHDLGIKYDP SIGIYGLDFYVVLGRPGFSI ADKKRRTGCIGAKHRISKEE AMRWFQQKYDGILPGK



Immunofluorescence analysis of MCF-7 cells using RPL11 antibody (STJ115503). Blue: DAPI for nuclear staining.