

Anti-NUTF2 antibody (1-127) (STJ115537)

STJ115537

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10204
Gene Symbol	NUTF2
Uniprot ID	NTF2_HUMAN
Immunogen	
Immunogen Region	1-127
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-127 of human NUTF2 (NP_005787.1).
Immunogen Sequence	MGDKPIWEQIGSSFIQHYYQ LFDNDRTQLGAIYIDASCLT WEGQQFQGKAAIVEKLSSLP FQKIQHSITAQDHPQTPDSC IISMVVGQLKADEDPIMGFH QMFLKNINDAWVCTNDMFR LALHNFG



Immunofluorescence analysis of U2OS cells using NUTF2 antibody (STJ115537). Blue: DAPI for nuclear staining.