

Anti-RAD54L2 antibody (1178-1467) (STJ115546)

STJ115546

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IHC-P:1:50-1:200
Range	IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23132
Gene Symbol	RAD54L2
Uniprot ID	ARIP4_HUMAN
Immunogen	
Immunogen Region	1178-1467
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1178-1467 of human RAD54L2 (NP_055921.2).
Immunogen Sequence	ELQQYADVAAARESRQSSPS TNAALPGPPAQLMDSSAVPG TALGTEPRLGGHCLNSSLLV TGQPCGDRHPVLDLRGHRK LATPPAAQESSRRRSRKGLH PAPVQPYEHGYPVSGGFAMP PVSLNHNLTTPFTSQAGENS LFMGSTPSYYQLSNLLADAR LVFPVTTDPLVPAGPVSSSS TATSVTASNPSFMLNPSVPG ILPSYSLPFSQPLLSEPRMF APFPSPVLPSNLSRGMSIY



Immunohistochemistry analysis of RAD54L2 in paraffin-embedded rat liver tissue using RAD54L2 rabbit polyclonal antibody (STJ115546) at a dilution of 1:200 (40x lens). High pressure antigen retrieval was performed with 0.01 M citrate buffer (pH 6.0) prior to immunohistochemistry staining.



Immunofluorescence analysis of MCF-7 cells using RAD54L2 antibody (STJ115546). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081