

Anti-SERPINA6 antibody (136-405) (STJ115569)
STJ115569

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	866
Gene Symbol	SERPINA6
Uniprot ID	CBG_HUMAN
Immunogen	
Immunogen Region	136-405
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 136-405 of human SERPINA6 (NP_001747.2).
Immunogen Sequence	SLELLESFSADIKHYEYEV LAMNFQDWATASRQINSYVK NKTQGKIVDLFSGLDSPAIL VLVNYIFFKGTWTQPFDLAS TREENFYVDETTVMKVPMMML QSSTISYLHDAELPCQLVQM NYVGNGTVFFILPDKGKMNT VIAALSRDTINRWSAGLTSS QVDLYIPKVTISGVYDLGDV LEEMGIADLFTNQANFSRIT QDAQLKSSKVVHKAVLQLNE EGVDTAGSTGVTLNLTSKP



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081