

Anti-EIF4G2 antibody (750-850) (STJ115575)

STJ115575

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1982
Gene Symbol	EIF4G2
Uniprot ID	IF4G2_HUMAN
Immunogen	
Immunogen Region	750-850
Specificity	A synthetic peptide corresponding to a sequence within amino acids 750-850 of human EIF4G2/p97 (NP_001409.3).
Immunogen Sequence	IYKWIKDNISPKLHVDKGFV NILMTSFLQYISSEVNPPSD ETDSSAPSKEQLEQEKQLL LSFKPVMQKFLHDHVDLQVS ALYALQVHCYNSNFPKGMLL R