

Anti-PROCR antibody (19-209) (STJ115600)

STJ115600

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10544
Gene Symbol	PROCR
Uniprot ID	EPCR_HUMAN
Immunogen	
Immunogen Region	19-209
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 19-209 of human PROCR (NP_006395.2).
Immunogen Sequence	QDASDGLQRLHMLQISYFRD PYHVVYQGNASLGGLHVL EGPDTNTTIQLQLQEPES WARTQSGLSYLLQFHGLVR LVHQERTLAFPLTIRCFGLG ELPPEGSRAHVFEVAVNGS SFVSFRPERALWQADTQVTS GVTFTLQQLNAYNRTRYEL REFLEDTCVQYVQKHISAEN TKGSQTSRSYT



Immunohistochemistry analysis of paraffin-embedded rat kidney using PROCR antibody (STJ115600) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded rat liver using PROCR antibody (STJ115600) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Western blot analysis of U-87MG, using PROCR antibody (STJ115600) at 1:600 dilution. Secondary antibody: HRP-Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081