

Anti-CD300LF antibody (23-150) (STJ115772)

STJ115772

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	146722
Gene Symbol	CD300LF
Uniprot ID	CLM1_HUMAN
Immunogen	
Immunogen Region	23-150
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 23-150 of human CD300LF (NP_001276011.1).
Immunogen Sequence	TQITGPTTVNGLERGS LTVQ CVYRSGWETYLKWWCRGAIW RDCKILVKTSGSEQEVKRDR VSIKDNQKNRTFTVTMEDLM KTDADTYWCGIEKTGNDLGV TVQVTIDPAPVTQEETSSSP TLTGHHLD



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081