

Anti-KRI1 antibody (410-709) (STJ115820)
STJ115820

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	65095
Gene Symbol	KRI1
Uniprot ID	KRI1_HUMAN
Immunogen	
Immunogen Region	410-709
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 410-709 of human KRI1 (NP_075384.3). EYYGAVEEEKPQFEIEEGLE DDWNWDTWDGPEQEGDWSQQ ELHCEDPNFNMADADYDPSQP RKKKKREAPLTGKKKRKSPFA
Immunogen Sequence	AAVGQEKPVFEPGDKTFEY LDEYYRLDYEDIIDDLPCR F KYRTVVPDFGLSTEEILAA DDKELNRWC SLKKTCTMYRSE QEELRDKRAYSQKAQNSWKK RQVFKSLCREEAETPAEATG KPQRDEAGPQRQLPALDGSL MGPESPPAQEEEEAPVSPHK



Western blot analysis of extracts of Jurkat cells, using KRI1 antibody (STJ115820) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081