

Anti-EWSR1 antibody (1-300) (STJ115913)

STJ115913

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:100
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2130
Gene Symbol	EWSR1
Uniprot ID	EWS_HUMAN
Immunogen	
Immunogen Region	1-300
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human EWSR1 (NP_001156757.1).
Immunogen Sequence	MASTDYSTYSQAAQGGYSA YTAQPTQGYAQTTQAYGQQS YGTYGQPTDVSYTQAQTTAT YGQTAYATSYGQPPTGYTTP TAPQAYSQPVQGYGTGAYDT TTATVTTTQASYAAQSAYGT QPAYPAYGQQPAATAPTRPQ DGNKPTETSQPSSTGGYNQ PSLGYGQSNYSYPQVPGSYP MQPVTAPPSYPPTSYSSTQP TSYDQSSYSQQNTYGQPSY GQQSSYGQQSSYGQQPPTS



Immunofluorescence analysis of U2OS cells using EWSR1 antibody (STJ115913).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081