

Anti-LGALS8 antibody (20-317) (STJ115954)

STJ115954

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3964
Gene Symbol	LGALS8
Uniprot ID	LEG8_HUMAN
Immunogen	
Immunogen Region	20-317
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-317 of human Galectin 8/Galectin 8/LGALS8 (NP_963838.1).
Immunogen Sequence	VTIPDQLDPGTLVIRGHV PSDADRFQVDLQNGSSMKPR ADVAFHFNPRFKRAGCIVCN TLINEKWGREEITYDTPFKR EKSFEIVIMVLKDKFQVAVN GKHTLLYGHRIQPEKIDTLG IYGVNIHSIGFSFSSDLQS TQASSLELTEISRENVKPSG TPQLRLPFAARLNTMPGPGR TVVVKGEVNANAKSFNVDLL AGKSKDIALHLNPRNLKAF VRNSFLQESWGEEERNITS