

Anti-RNF8 antibody (1-260) (STJ116018)

STJ116018

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/CC:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9025
Gene Symbol	RNF8
Uniprot ID	RNF8_HUMAN
Immunogen	
Immunogen Region	1-260
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human RNF8 (NP_003949.1).
Immunogen Sequence	MGEPGFFVTGDRAGGRSWCL RRVGMSAGWLLLEDGCEVTV GRGFGVTYQLVSKICPLMIS RNHCVLKQNPGEQWTIMDNK SLNGVWLNRRARLEPLRVYSI HQGDYIQLGVPLENKENAEY EYEVTEEDWETIYPCLSPKN DQMIKKNKELRTKRKFSLDE LAGPGAEGPSNLKSKINKVS CESGQPVKSQKGGEVASTPS DNLDPKLTALPESKTTGAPI YPGFPKVTEVHHEQKASNS



Immunofluorescence analysis of MCF-7 cells using RNF8 antibody (STJ116018). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081